



Nidhi Meditech System Provides genuine endo urological disposables like Double Pigtail Ureteral stent, Mono 'J' Stent, Endopylotomy stent, Ureteric cathaters etc. All the endourological disposables manufactured & as per ISO 13485 : 2003 and CE 1023 standards.

All the products manufactured from tecoflex grade Polyurethanes according to ISO-10993 series

In general, the ease of processing and high strength of polyurethane makes them the material of choice in medical devices. Polyurethanes pass the tests of biocompatibility and chemical stability. They can be easily loaded with radiopaque material for X-ray detection. However, the major advantages of polyurethane are their thermoliability that makes them very soft inside the body, minimizing the sense of discomfort associated with implantable foreign body. Nidhi double loop stent is made of tecoflex, which has the maximum degree of softening inside the body. It has a wide margin of durometer extending from very soft to very hard with different sizes and length according to the physician specification.

DOUBLE PIGTAIL URETERAL STENT

FEATURES

- 1. Durable with superior bio-compatibility
- 2. Thin wall (provides potency and good drainage).
- 3. Astraumatic insertion is revealed quite smooth.
- 4. Minimized tissue irritation with easy but rapid astraumatic insertion.
- 5. Cost effective.

DOUBLE PIGTAIL URETERAL STENT WITH PUSHER (WHITE / BLUE COLOR)

Size	Product No.		Length	Compatible Guide Wire
4.5 FR	NMS-DJ-4.5 FR	26 (CE)	26 cms.	0.025"
4.5 FR	NMS-DJ-4.5 FR	26 (BEO)	26 cms.	0.025"
4.8 FR	NMS-DJ-4.8 FR	26 (CE)	26 cms.	0.032"
4.8 FR	NMS-DJ-4.8 FR	26 (BEO)	26 cms.	0.032"
5 FR	NMS-DJ- 5 FR	26 (CE)	26 cms.	0.035"
5 FR	NMS-DJ-5 FR	26 (BEO)	26 cms.	0.035"
5.5 FR	NMS-DJ- 5.5 FR	26 (CE)	26 cms.	0.035"
5.5 FR	NMS-DJ-5.5 FR	26 (BEO)	26 cms.	0.035"
6 FR	NMS-DJ- 6 FR	26 (CE)	26 cms.	0.038"
6 FR	NMS-DJ- 6 FR	26 (BEO)	26 cms.	0.038"

DOUBLE PIGTAIL URETERAL STENT WITH PUSHER & SUTURE (WHITE / BLUE COLOR)


	,	,		
Size	Product No.		Length	Compatible Guide Wire
4.5 FR	NMS-DJS-4.5 FR	26 (CE)	26 cms.	0.025"
4.5 FR	NMS-DJS-4.5 FR	26 (BEO)	26 cms.	0.025"
4.8 FR	NMS-DJS-4.8 FR	26 (CE)	26 cms.	0.032"
4.8 FR	NMS-DJS-4.8 FR	26 (BEO)	26 cms.	0.032"
5 FR	NMS-DJS- 5 FR	26 (CE)	26 cms.	0.035"
5 FR	NMS-DJS-5 FR	26 (BEO)	26 cms.	0.035"
5.5 FR	NMS-DJS- 5.5 FR	26 (CE)	26 cms.	0.035"
5.5 FR	NMS-DJS-5.5 FR	26 (BEO)	26 cms.	0.035"
6 FR	NMS-DJS- 6 FR	26 (CE)	26 cms.	0.038"
6 FR	NMS-DJS- 6 FR	26 (BEO)	26 cms.	0.038"

Please specify in your order color = White/Blue, product code, Single Loop / Multi Loop

Any special size & length can be available on request.


Application	Applied for internal drainage from the uretropelvic
	junction to the bladder, as a provisional measure.
Features	Made of soft telecoflex polyuretharie for
	Smooth Passage radio-opeque implant.
Use	For single use
Packing	Tyvek pouch as per EN-868-5-series standard
Sterilizations	With Ethylene Oxide EN-550 series standard

CE = One End Close, BEO = Both End Open





DOUBLE PIGTAIL URETERAL STENT SET



Application	Applied for internal drainage from the uretropelvic
	junction to the bladder, as a provisional measure.
Features	Made of soft telecoflex polyuretharie for
	Smooth Passage radio-opeque implant.
Contents	Pusher, Suture, Clamp with PTFE coated Guidewire
Use	For single use
Packing	Tyvek pouch as per EN-868-5-series standard
Sterilizations	With Ethylene Oxide EN-550 series standard
Packing Area	Class 100 as per ISO 14644-1

Size	Product No.		Length	Compatible Guide Wire
4.5 FR	NMS-DJ-ST-4.5 FR	26 (CE)	26 cms.	0.025"
4.5 FR	NMS-DJ-ST-4.5 FR	26 (BEO)	26 cms.	0.025"
4.8 FR	NMS-DJ-ST-4.8 FR	26 (CE)	26 cms.	0.032"
4.8 FR	NMS-DJ-ST-4.8 FR	26 (BEO)	26 cms.	0.032"
5 FR	NMS-DJ-ST- 5 FR	26 (CE)	26 cms.	0.035"
5 FR	NMS-DJ-ST-5 FR	26 (BEO)	26 cms.	0.035"
5.5 FR	NMS-DJ-ST- 5.5 FR	26 (CE)	26 cms.	0.035"
5.5 FR	NMS-DJ-ST-5.5 FR	26 (BEO)	26 cms.	0.035"
6 FR	NMS-DJ-ST- 6 FR	26 (CE)	26 cms.	0.038"
6 FR	NMS-DJ-ST- 6 FR	26 (BEO)	26 cms.	0.038"

URETERIC CATHATER



C€ 1023

CE = One End Close, BEO = Both End Open

URETERIC CATHETER (OPEN END)

Application	Used for drainage as well as retrograde pyelogram
	by open end tip it is placed over a guide wire for
	negotiation to tortuous ureter.
Features	Smooth passage in tortuous ureter.
Use	For single use
Packing	Tyvek pouch as per EN-868-5-series standard
Sterilizations	With Ethylene Oxide EN-550 series standard
Packing Area	Class 100 as per ISO 14644-1

Siz	e Pro	duct No.		Length	
3 F	R NM	S-UC-3 FR	70 (OE)	70 cms.	
4 F	R NM	S-UC-4 FR	70(OE)	70 cms.	
5 F	R NM	S-UC-5 FR	70 (OE)	70 cms.	
6 F	R NM	S-UC-6 FR	70(0E)	70 cms.	

URETERIC CATHETER (CLOSED END)

OHETEHIO OATHETEH (OLOGED END)			
Application	Used for drainage and retrograde pyelogram.		
Features	Best in its inherent quality providing drainage and		
	retrograde pyelogram.		
Use For single use Packing Tyvek pouch as per EN-868-5-series standard Sterilizations With Ethylene Oxide EN-550 series standard			
		Packing Area	Class 100 as per ISO 14644-1

Size	Product No.	Length
3 FR	NMS-UC-3 FR 70 (CE)	70 cms.
4 FR	NMS-UC-4 FR 70 (CE)	70 cms.
5 FR	NMS-UC-5 FR 70 (CE)	70 cms.
6 FR	NMS-UC-6 FR 70 (CE)	70 cms.



MONO-J-STENT





Application By direct puncture, it provides temporary drainage of

urine, from the kidney. This is done in presence of

distal obstruction.

Features Provides drainage from kidney, smoothly.

Use For single use

Packing Tyvek pouch as per EN-868-5-series standard
Sterilizations With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

MONO-J-STENT

Size	Product No.	Length
6 FR	NMS-PCNC-6 FR- 22	22 cms.
6 FR	NMS-PCNC-6 FR 30	30 cms.
8 FR	NMS-PCNC-8 FR 22	22 cms.
8 FR	NMS-PCNC-8 FR - 30	30 cms.
9 FR	NMS-PCNC-9 FR 22	22 cms.
9 FR	NMS-PCNC-9 FR 30	30 cms.
10 FR	NMS-PCNC-10FR 22	22 cms.
10 FR	NMS-PCNC-10FR 30	30 cms.
12 FR	NMS-PCNC-12FR 22	22 cms.
12 FR	NMS-PCNC-12FR 30	30 cms.
14 FR	NMS-PCNC-14FR 22	22 cms.
14 FR	NMS-PCNC-14FR 30	30 cms.

MONO-J-STENT WITH TROCAR

Size	Product No.	Length		
6 FR	NMS-PCNC-N-6 FR- 22	22 cms.		
6 FR	NMS-PCNC-N-6 FR 30	30 cms.		
8 FR	NMS-PCNC-N-8 FR 22	22 cms.		
8 FR	NMS-PCNC-N-8 FR - 30	30 cms.		
9 FR	NMS-PCNC-N-9 FR 22	22 cms.		
9 FR	NMS-PCNC-N-9 FR 30	30 cms.		
10 FR	NMS-PCNC-N-10FR 22	22 cms.		
10 FR	NMS-PCNC-N-10FR 30	30 cms.		
12 FR	NMS-PCNC-N-12FR 22	22 cms.		
12 FR	NMS-PCNC-N-12FR 30	30 cms.		
14 FR	NMS-PCNC-N-14FR 22	22 cms.		
14 FR	NMS-PCNC-N-14FR 30	30 cms.		

MONO-J-STENT KIT





MONO-J-STENT SET

C€ 1023



By direct puncture, it provides temporary drainage of

urine, from the kidney. This is done in presence of

distal obstruction.

Features Provides drainage from kidney, smoothly.

Contents Mono-J-Stent, Trocar, Guide Wire, Clamp, Scalpel

Use For single use

Packing Tyvek pouch as per EN-868-5-series standard Sterilizations With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

ENDOPYLOTOMY STENT

ENDOPYLOTOMY STENT

Application Used for Stenting the ureteropelvic junction & ureter

following a percutaneous endopyelotomy while providing both nephrostomy & ureteral drainage

Use For single use

PackingTyvek pouch as per EN-868-5-series standardSterilizationsWith Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

C€ 1023

5 fr / 10 fr 26 cms NMS-EPS-5/10FR26

6 fr / 12 fr 26 cms NMS-EPS-6/12FR26

7 fr / 14 fr 26 cms NMS-EPS-7/14FR26



MALECOT NEPHROSTOMY CATHETER

Malecot Nephrostomy Catheter

10 fr	30 cms	NMS-MC-10FR-30
12 fr	30 cms	NMS-MC-12FR-30
14 fr	30 cms	NMS-MC-14FR-30

Malecot Nephrostomy Catheter (With Needle)

10 fr	30 cms	NMS-MC-N-10FR-30
12 fr	30 cms	NMS-MC-N-12FR-30
14 fr	30 cms	NMS-MC-N-14FR-30





Application The Malecot nephrostomy catheter is placed percutaneously through an existing

nephrostomy tract. It is used for nephrostomy drainage in the renal pelvis.

Features Smooth passage in tortuous ureter.

Use For single use

Packing Tyvek pouch as per EN-868-5-series standard Sterilizations With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

STONE RETRIEVAL BASKET

Stone Retrieval Basket - Helical (S.S.)



Application To remove multiple stone fragments from ureter with

semi rigid endoscope.

Features Maximum patient safety and compatibility. easy to operate.

Use For single use

Packing Blister tray with tyvek paper as per EN-868-5-series standard

Sterilizations With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

Size (In FR)	Product No.	Length
3 / 4 WIRE	NMS-DB-H-3FR-4FR-70	70 cms.
4 / 4 WIRE	NMS-DB-H-4FR-4FR-70	70 cms.
4 / 5 WIRE	NMS-DB-H-4FR-5FR-70	70 cms.
5 / 4 WIRE	NMS-DB-H-5FR-4FR-70	70 cms.
5 / 5 WIRE	NMS-DB-H-5FR-5FR-70	70 cms

Flatwire Stone Retrieval Basket - Segura (S.S.)



Application It is used for removing the large

stones from ureter with semi rigid endoscope.

Features Retrieving the large stones, smoothly.

Use For single use

Packing Blister tray with tyvek paper as per EN-868-5-series standard

Sterilizations With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

Size (In FR)	Product No.	Length
3 / 4 WIRE	NMS-DB-S-3FR-4FR-70	70 cms.
4 / 4 WIRE	NMS-DB-S-4FR-4FR-70	70 cms.
4 / 5 WIRE	NMS-DB-S-4FR-5FR-70	70 cms.
5 / 4 WIRE	NMS-DB-S-5FR-4FR-70	70 cms.
5 / 5 WIRE	NMS-DB-S-5FR-5FR-70	70 cms.



STONE RETRIEVAL BASKET

Flatwire Stone Retrieval Basket - '0' Tip (S.S.)



Application To remove multiple stone fragments from ureter with

semi rigid endoscope.

Features Maximum patient safety and compatibility. Easy to operate.

Use For single use

Packing Blister tray with tyvek paper as per EN-868-5-series standard

Sterilizations With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

Size (In FR)	Product No.	Length
3 / 4 WIRE	NMS-DB-H-3FR-4FR-70	70 cms.
4 / 4 WIRE	NMS-DB-H-4FR-4FR-70	70 cms.
4 / 5 WIRE	NMS-DB-H-4FR-5FR-70	70 cms.
5 / 4 WIRE	NMS-DB-H-5FR-4FR-70	70 cms.
5 / 5 WIRE	NMS-DB-H-5FR-5FR-70	70 cms.

Flatwire Stone Retrieval Basket - Segura (Nitinol)



Application It is used for removing the large stones from Ureter with

semi rigid/flexible endoscope.

Features Retrieving the large stones, smoothly.

Use For single use

Packing Blister tray with tyvek paper as per EN-868-5-series standard

Sterilizations With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

Size (In FR)	Product No.	Length
3 / 4 WIRE	NMS-DB-S-3FR-4FR-70	70 cms.
4 / 4 WIRE	NMS-DB-S-4FR-4FR-70	70 cms.
4 / 5 WIRE	NMS-DB-S-4FR-5FR-70	70 cms.
5 / 4 WIRE	NMS-DB-S-5FR-4FR-70	70 cms.
5 / 5 WIRE	NMS-DB-S-5FR-5FR-70	70 cms.

Stone Retrieval Basket - Helical (Nitinol)



Application To remove multiple stone fragments from ureter with

semi rigid/flexible endoscope.

Features Maximum patient safety and compatibility. easy to operate.

Use For single use


Packing Blister tray with tyvek paper as per EN-868-5-series standard

Sterilizations With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

Size (In FR)	Product No.	Length
3 / 4 WIRE	NMS-DB-H-3FR-4FR-70	70 cms.
4 / 4 WIRE	NMS-DB-H-4FR-4FR-70	70 cms.
4 / 5 WIRE	NMS-DB-H-4FR-5FR-70	70 cms.
5 / 4 WIRE	NMS-DB-H-5FR-4FR-70	70 cms.
5 / 5 WIRE	NMS-DB-H-5FR-5FR-70	70 cms.

Flatwire Stone Retrieval Basket - Segura (Nitinol)



Application It is used for removing the large

stones from ureter with semi rigid/flexible endoscope.

Features Retrieving the large stones, smoothly.

Use For single use

Packing Blister tray with tyvek paper as per EN-868-5-series standard

Sterilizations With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

Size (In FR)	Product No.	Length
3 / 4 WIRE	NMS-DB-S-3FR-4FR-70	70 cms.
4 / 4 WIRE	NMS-DB-S-4FR-4FR-70	70 cms.
4 / 5 WIRE	NMS-DB-S-4FR-5FR-70	70 cms.
5 / 4 WIRE	NMS-DB-S-5FR-4FR-70	70 cms.
5 / 5 WIRE	NMS-DB-S-5FR-5FR-70	70 cms.



SINGLE STEP NOTTINGHAM DIALATOR

Aplication: Used to dilate the intramural tunnel and ureter in one step.

May also be used to inject contrast medium to determine sub mucosal tunneling or anesthetic jelly to aid in dislodging

ureteral stones

Features - Smooth surfaces for ease of introduction.

> - Uniform taper to reduce trauma. - Used with 0.038" guide wire

Size 10 fr, 12 fr & 14 fr

Length 70 cm Use For single use

Packing Tyvek pouch as per EN-868-5-series standard **Sterilizations** With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

FASCIAL DIALATORS SET (PCN)



Aplication: Used for fasecial dilation of the ureter prior to

ureteroscopy and / or stone manipulation

Features - Smooth surfaces for ease of introduction.

- Uniform taper to reduce trauma.

Size 6 fr, 16 fr Length Use

For single use

Packing Tyvek pouch as per EN-868-5-series standard **Sterilizations** With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1

URETERAL DIALATORS SET



Used for dilation of ureter prior to uretero renoscopy and

also for sequential (step by step) dilation of the intramural tunnel and ureter and / or stone manipulation.

Features

- Smooth surfaces for ease of introduction.

- Uniform taper to reduce trauma.

Size 16 fr Length 70 cm Use For single use

Packing Tyvek pouch as per EN-868-5-series standard Sterilizations With Ethylene Oxide EN-550 series standard

Packing Area Class 100 as per ISO 14644-1





(Complete Urological Solutions Since-1995)

Regd. Off Add:

A/207, Tirthraj Complex, Opp. V.S.Hospital, Ellisbridge, Ahmedabad 380006 (India)

Tel: ++91-79-26578422, 26578322 Fax: ++91-79-30028422

Cell: +91-9825070980 / 9328270980

Email: info@nidhimeditech.com, sales@nidhimeditech.com

Website: www.nidhimeditech.com

www.nidhimeditech.com